

# Protein Quantitation II: Multiple Reaction Monitoring


Kelly Ruggles

[kelly@fenyolab.org](mailto:kelly@fenyolab.org)


New York University

# Traditional Affinity-based proteomics


## Use antibodies to quantify proteins




Western Blot



ELISA




Immunohistochemistry




RPPA



Immunofluorescence




# Mass Spectrometry based proteomic quantitation




# Multiple Reaction Monitoring (MRM)

- Triple Quadrupole acts as ion filters
- Precursor selected in first mass analyzer (Q1)
- Fragmented by collision activated dissociation (Q2)
- One or several of the fragments are specifically measured in the second mass analyzer (Q3)




# Peptide Identification with MRM




- Transition: Precursor-Fragment ion pair are used for protein identification
- Select both Q1 and Q3 prior to run
  - Pick Q3 fragment ions based on discovery experiments, spectral libraries
  - Q1 doubly or triply charged peptides
- Use the 3 most intense transitions for quantitation

# Label-free quantification

- Usually use 3 or more precursor-product ion pairs (transitions) for quantitation
- Relies on direct evaluation of MS signal intensities of naturally occurring peptides in a sample.
- Simple and straightforward
- Low precision
- Several peptides for each protein should be quantified to avoid false quantification



# Stable Isotope Dilution (SID)



- Use isotopically labeled reference protein
- $^{13}\text{C}$  and/or  $^{15}\text{N}$  labeled peptide analogs
- Chemically identical to the target peptide but with mass difference
- Add known quantity of heavy standard
- Compare signals for the light to the heavy reference to determine for precise quantification

Q1                      Q3




## Heavy





# Quantification Details



SIS: Stable Isotope Standard

PAR: Peak Area Ratio

$$\text{PAR} = \frac{\text{Light (Analyte) Peak Area}}{\text{Heavy (SIS) Peak Area}}$$

Analyte concentration = PAR \* SIS peptide concentration

- Use at least 3 transitions
- Have to make sure these transitions do not have interferences

# Strengths of MRM

- Can detect multiple transitions on the order of 10msec per transition
- Can analyze many peptides (100s) per assay and the monitoring of many transitions per peptide
- High sensitivity
- High reproducibility
- Detects low level analytes even in complex matrix
- Golden standard for quantitation!


# Weaknesses of MRM

- Focuses on defined set of peptide candidates
  - Need to know charge state, retention time and relative product ion intensities before experimentation
- Physical limit to the number of transitions that can be measured at once
  - Can get around this by using time-scheduled MRM, monitor transitions for a peptide in small window near retention time


# Parallel Reaction Monitoring (PRM)

- Q3 is substituted with a high resolution mass analyzer to detect all target product ions
- Generates high resolution, full scan MS/MS data
- All transitions can be used to confirm peptide ID
- Don't have to choose ions beforehand

## A SRM




## B PRM



# SWATH-MS: Data Collection


- Data acquired on quadrupole-quadrupole TOF high resolution instrument cycling through 32-consecutive 25-Da precursor isolation windows (swaths).
- Generates fragment ion spectra for all precursor ions within a user defined precursor retention time and  $m/z$
- Records the fragment ion spectra as complex fragment ion maps




32 discrete precursor isolation windows of 25-Da width across the 400-1200  $m/z$  range

# Applications of MRM


Metabolic pathway analysis




Protein complex  
subunit stoichiometry




Phosphorylation



Modifications within protein




Biomarkers: protein indicator  
correlating to a disease state



# MRM and Biomarker Verification


- Measurable indicator that provides the status of a biological state
  - Diagnosis
  - Prognosis
  - Treatment efficacy
- Shotgun proteomics → Biomarker Discovery (<100 patients)
- Targeted proteomics → Biomarker Validation (~1000s patients)
  - Requires higher threshold of certainty
  - Remove high false positives from discovery phase
- Most often plasma/serum, but can be tissue-based biomarkers




Meng Z and Veenstra TD, 2011

# MRM and Biomarker Verification

- Originally used to analyze small molecules since the late 1970s
- More recently, used for proteins and peptide quantitation in complex biological matrices
- With small molecules, the matrix and analyte have different chemical natures so separation step is able to remove other components from analytes




- With proteomics, both the analytes and the background matrix are made up of peptides, so this separation cannot occur. Leads to decreased sensitivity and increased interference.





# Enhancing MRM Sensitivity for Biomarker Discovery


## Sample Enrichment



Shi T., et al. 2012

## MRM3


Further fragments product ions  
Reduces background




Meng Z and Veenstra TD, 2011

# Workflow of MRM and PRM MS/MS

## SRM




## PRM




# Workflow of MRM and PRM MS/MS

## SRM




## PRM



Define a set of proteins based on clinical/biological question


# Motivating Example: AKT1 and Breast Cancer




- AKT
- PDK
- BAD
- MDM2
- GSK3
- mTOR
- RAF1

# Workflow of MRM and PRM MS/MS

## SRM



## PRM




- Proteotypic
- Consistently observed by LC-MS methods

# Selecting Peptides


- A few representative peptides will be used to quantify each protein
- Need to fulfill certain characteristics
  - Have an unique sequence
  - Consistently observed by LC-MS methods
  - 8-25 amino acids
  - Good ionization efficiency
  - $m/z$  within the range of the instrument
  - No missed cleavages
  - Not too hydrophilic (poorly retained) or hydrophobic (may stick to column)

# Identifying Proteotypic Peptides



# LC/MS Properties: GPMDB


- Compares peptides to a collection of previously observed results
- Determines how many times the peptide has been observed by others
- Most proteins show very reproducible peptide patterns






# LC/MS Properties: Skyline

- Compares peptides to MS/MS spectral library
- Predicts most abundant transitions




# Workflow of MRM and PRM MS/MS

## SRM




## PRM



PRM allows for selection of transitions post-data acquisition


# Selecting Transitions

- Limitation of MRM-MS:  $\sim 1$ -2  $m/z$  unit window for precursor and fragment ion occasionally let in interfering peptides with similar characteristics
- If we want to use these transitions for quantitation, we need to be confident there are no interferences
- Largest always largest, smallest always smallest etc.
- b-fragments of high  $m/z$  are less represented on QqQ



# Selecting Transitions

- Limitation of MRM-MS:  $\sim 1$ -2  $m/z$  unit window for precursor and fragment ion occasionally let in interfering peptides with similar characteristics
- If we want to use these transitions for quantitation, we need to be confident there are no interferences
- Largest always largest, smallest always smallest etc.
- b-fragments of high  $m/z$  are less represented on QqQ



# Selecting Transitions: SRMCollider

- Input peptides of interest
- Determines the m/z values for transition pair
- Simulates a typical SRM experiment
- Predicts fragment intensities and retention time information for input peptide
- Compares the transition to all other transitions in a background proteome
- Outputs the number of predicted interferences for each transition for that peptide

SRM Collider

version 1.4  
Hannes Röst 2012

Collider

Download

About

Instructions

The SRM Collider is a program that will take your input transitions and compare them to all other transitions in a given background proteome and find interferences. It will report these interferences on a per-peptide basis, allowing a researcher to identify peptides that share many transitions with the target peptide.

Please enter the peptide sequences here (see [Instructions](#) for help):

YDEGMDCMDNER

**Input peptide sequence**

SSRCalc window

Q1 mass window

Q3 mass window

Low mass threshold for transitions

High mass threshold for transitions

Genome

Consider isotopes up to

Missed Cleavages

Find ULS up to order\*

Charge check: ☒ Check that interfering signals

Modifications: ☐ oxidized Methionines ☐

**Peptide YDEGMDCMDNER**

Sequence	Q1	Q1 window	SSRCalc	SSRCalc window	Interfering precursors	Background	Graph
YDEGMDCMDNER	796.769077374	± 0.35	15.82	± 1.0	32	human	<a href="#">Graph</a>



**Transition Overview**

Transition	Q3	Interferences	Graph
y10	1185.4	0	<a href="#">Graph</a>
y9	1070.37	0	<a href="#">Graph</a>
y6	767.28	0	<a href="#">Graph</a>
y5	664.27	0	<a href="#">Graph</a>
b3	408.14	0	<a href="#">Graph</a>
b4	523.17	0	<a href="#">Graph</a>
b5	580.19	0	<a href="#">Graph</a>
b6	711.23	0	<a href="#">Graph</a>
b7	826.26	0	<a href="#">Graph</a>
b10	1175.33	0	<a href="#">Graph</a>
y11	1314.44	1	<a href="#">Graph</a>
y8	1013.35	1	<a href="#">Graph</a>
y7	882.31	1	<a href="#">Graph</a>
y4	533.23	1	<a href="#">Graph</a>
y3	418.21	1	<a href="#">Graph</a>
y2	304.16	1	<a href="#">Graph</a>
b8	929.27	1	<a href="#">Graph</a>
b9	1060.31	1	<a href="#">Graph</a>
y12	1429.47	2	<a href="#">Graph</a>
b11	1289.38	3	<a href="#">Graph</a>
b12	1418.42	3	<a href="#">Graph</a>

**Choose peptides that have at least one transition with zero interferences**

# Selecting Transitions: Skyline


- Can use to find best transitions to pick
  - Intensity (rank)
  - Dot product (similarity to reference spectra)




Want high rank and dotp close to 1

# Workflow of MRM and PRM MS/MS

## SRM





## PRM




# Validating Transitions: Contrast Angle

- Spectral Contrast Angle: each spectrum represented as a vector in N-dimensional space
- Spectra that resemble each other have vectors pointing in the same direction ( $\theta \sim 0^\circ$ )



$$\cos\theta = \frac{\sum a_i b_i}{\sqrt{\sum a_i^2 \cdot \sum b_i^2}}$$



$$r_a = \sqrt{a_i^2}$$

$$r_b = \sqrt{b_i^2}$$




# Validating Transitions: “Branching ratio”

Branching Ratio (BR): ratio of the peak intensities


$$BR = \ln \left\{ \frac{\frac{I_{Ax}}{I_{Bx}}}{\frac{\sum \frac{I_{AxS}}{I_{BxS}}}{n}} \right\}$$

$I_{Ax}$ ,  $I_{Bx}$  : Peak areas of Analyte  
 $I_{AxS}$ ,  $I_{BxS}$  : Peak areas of SIS  
 $n$ =number of SIS transitions


Light (Analyte)



Heavy(SIS)




# Validating Transitions in MRM: AuDIT



- AuDIT: Automated Detection of Inaccurate and imprecise Transitions
- Uses “branching ratio”
  1. Calculate relative ratios of each transition from the same precursor
  2. Apply t-test to determine if relative ratios of analyte are different from relative ratios of SIS

# Validating Transitions in MRM: AuDIT




Blue: Light  
Red: Heavy

Relative product ions should have a constant relationship

# Validating Transitions in PRM: CRAFTS



- PRM and MRM are most useful when quantifying protein in a complex matrix
  - Tumor lysate
  - Plasma
- Simple Matrix (buffer) should have no interferences
- Compare the transitions in complex to those in simple
- Ratio close to 1 indicates low interference

# Validating Transitions in PRM: CRAFTS



- Simple matrix: peptide carrier solution
- Complex matrix: unfractionated tumor digest
- Simple matrix should have minimal interference- use this as reference
- Transitions in complex buffer should have the same relative intensities of transitions within the spectra
- Transitions in complex with relative intensities different from simple → interference

# Validating Transitions in PRM: CRAFTS



Transition	Simple	Complex
y2	100	150
y5	40	220
y10	60	90

Ratio of  
Transitions

	y2	y5	y10
y2	1	$y5/y2$	$y10/y2$
y5	$y2/y5$	1	$y10/y5$
y10	$y2/y10$	$y5/y10$	1

Simple  
Matrix

	y2	y5	y10
y2	1	0.4	0.6
y5	2.5	1	1.5
y10	1.67	0.67	1

Complex  
Matrix

	y2	y5	y10
y2	1	<b>1.47</b>	0.6
y5	<b>0.68</b>	1	<b>0.41</b>
y10	1.67	<b>2.44</b>	1

# Validating Transitions in PRM: CRAFTS

Ratio of  
Transitions

	y2	y5	y10
y2	1	y5/y2	y10/y2
y5	y2/y5	1	y10/y5
y10	y2/y10	y5/y10	1

Simple  
Matrix

	y2	y5	y10
y2	1	<b>0.4</b>	0.6
y5	<b>2.5</b>	1	<b>1.5</b>
y10	1.67	<b>0.67</b>	1


Complex  
Matrix

	y2	y5	y10
y2	1	<b>1.47</b>	0.6
y5	<b>0.68</b>	1	<b>0.41</b>
y10	1.67	<b>2.44</b>	1

Complex/Simple

	y2	y5	y10
y2	1	<b>3.675</b>	1
y5	<b>0.272</b>	1	<b>0.273</b>
y10	1	<b>3.641</b>	1


# Validating Transitions in PRM: CRAFTS






# Validating Transitions in PRM: CRAFTS


Light




Heavy



Minimal Interference: y3, y4, b3, y5, y6, y7, b6, b7, y8  
With Interference: y9, y11




Minimal Interference: y3, y4, b3, y5, y6, y7, y10, y11  
With Interference: y8, y9




Use highlighted values to get mean ratio


# CRAFTS: Ranking Transitions by Mean Combinatorial Ratio




## Peak Area Ratios in Simple and Complex Matrix




## Branching Ratio Matrices




## Combinatorial Ratios




## Normalized Combinatorial Ratios




## Apply Threshold (T1)




## Identify Transitions without Interference



## Rank by Mean Combinatorial Ratio




# Open Source MRM analysis tools




# SKYLINE for creating targeted MS/MS methods

Skyline digests proteins and fragments peptides and uses spectral library to find transition intensity



# Skyline for MRM: Method Building



Input all peptides of interest


Shows graphs of MS/MS spectra from spectral library

# Skyline for MRM: Method Building

- Helps generate prototypic peptide lists using MS/MS spectral libraries
- Find which peptides can be measured in specific matrix
- Find best transitions to measure for a peptide
- Creates transition lists and vendor-specific instrument methods for MRM experiments

# Skyline for MRM: Quantification


- Import raw files into skyline
- Pick peptide of interest
- Check standard peaks





# Skyline for MRM: Quantification


- Use the heavy standard PAR to make calibration curve
- Determine sample quantity based on curve



Questions?

# MRM Instrumentation

## Triple Quadrupole



## Quadrupole Time-of-Flight (Qqtof)

