


Protein quantitation I: Overview (Week 6)


Proteomic Bioinformatics - Quantitation



$$I_{ik} = \alpha_k \sum_j \left[C_{ij} p_{ij}^L p_{ij}^{Pr} p_{ijk}^D \right] p_{ik}^{Pep} p_{ik}^{LC} p_{ik}^{MS}$$

$$C_{ij}^k = \frac{I_{ik}}{\alpha_k p_{ij}^L p_{ij}^{Pr} p_{ijk}^D p_{ik}^{Pep} p_{ik}^{LC} p_{ik}^{MS}}$$

Quantitation - Label-Free (Standard Curve)



Quantitation - Label-Free (MS)


Sample i
Protein j
Peptide k

Assumption:


$$\alpha_k p_{ij}^L p_{ij}^{\text{Pr}} p_{ijk}^D p_{ik}^{\text{Pep}} p_{ik}^{\text{LC}} p_{ik}^{\text{MS}}$$

constant for all samples


$$C_{i_n j} / C_{i_m j} = I_{i_n j} / I_{i_m j}$$



Quantitation - Metabolic Labeling



Comparison of metabolic labeling and label-free quantitation



Label free assumption:

$$\alpha_k p_{ij}^L p_{ij}^{\text{Pr}} p_{ijk}^D p_{ik}^{\text{Pep}} p_{ik}^{\text{LC}} p_{ik}^{\text{MS}}$$


constant for all samples

Metabolic labeling assumption:

$$p_{ij}^M$$

constant for all samples and
the behavior of heavy and light
isotopes is identical

Intensity variation between runs




Replicates


1 IP
1 Fractionation
1 Digestion

VS

3 IP
3 Fractionations
1 Digestion




How significant is a measured change in amount?



It depends on the size of the random variation of the amount measurement that can be obtained by repeat measurement of identical samples.


Protein Complexes




↓
Digestion

↓
Mass spectrometry

Protein Complexes - specific/non-specific binding



Protein Turnover




$$\frac{dC_j^H(t)}{dt} = -(K_c + K_T)C_j^H(t) \quad \Rightarrow \quad C_j^H(t) = C_j^H(0)e^{-(K_c + K_T)t}$$

$$C_j^L(t) + C_j^H(t) = C_j^H(0)$$


$K_c = \log(2)/t_c$, t_c is the average time it takes for cells to go through the cell cycle, and $K_T = \log(2)/t_T$, t_T is the time it takes for half the proteins to turn over.

$$\log\left(\frac{I_j^H(t) + I_j^L(t)}{I_j^H(t)}\right) = t\left(\frac{1}{t_c} + \frac{1}{t_T}\right)\log(2)$$

Super-SILAC




Quantitation - Protein Labeling




Assumption: All losses after mixing are identical for the heavy and light isotopes and

$$p_{i_n j}^L p_{i_n j}^M \approx p_{i_m j}^L p_{i_m j}^M$$




Quantitation - Labeled Proteins




Assumption: All losses after mixing are identical for the heavy and light isotopes and


$$p_{i_n j}^L p_{i_n j}^M \approx p_{i_m j}^M$$



Quantitation - Labeled Chimeric Proteins




Quantitation - Peptide Labeling




Assumption: All losses after mixing are identical for the heavy and light isotopes and

$$p_{i_n j}^L p_{i_n j}^{\text{Pr}} p_{i_n jk}^D p_{i_n k}^M \approx$$

$$\approx p_{i_m j}^L p_{i_m j}^{\text{Pr}} p_{i_m jk}^D p_{i_m k}^M$$



Quantitation - Labeled Synthetic Peptides




Enrichment with
Peptide antibody

Anderson, N.L., et al.
Proteomics 3 (2004) 235-44

Assumption: All losses after mixing are identical for the heavy and light isotopes and


$$p_{i_n^j}^L p_{i_n^j}^{\text{Pr}} p_{i_n^{jk}}^D p_{i_n^k}^M \approx p_{sk}^M$$

**Synthetic
Peptides
(Heavy)**




Gerber et al. PNAS 100 (2003) 6940



Quantitation - Label-Free (MS/MS)




Quantitation - Labeled Synthetic Peptides




Quantitation - Isobaric Peptide Labeling




Quantitation – Label-Free (MS)




Quantitation – Label-Free (MS/MS)




Quantitation – Label-Free (Standard Curve)




Quantitation – Metabolic Labeling




Quantitation – Protein Labeling




Quantitation – Labeled Chimeric Proteins




Quantitation – Peptide Labeling




Quantitation – Isobaric Peptide Labeling







Quantitation – Labeled Synthetic Peptides




Isotope distributions



Isotope distributions




Estimating peptide quantity







Time dimension

Time

Intensity



m/z




Time





m/z

Sampling




Sampling




Acquisition time = 0.05σ



Sampling



Retention Time Alignment



Estimating peptide quantity by spectrum counting



What is the best way to estimate quantity?

Peak height

- resistant to interference
- poor statistics

Peak area

- better statistics
- more sensitive to interference





Curve fitting


- better statistics
- needs to know the peak shape
- slow

Spectrum counting


- resistant to interference
- easy to implement
- poor statistics for low-abundance proteins

Examples - qTOF








Examples - Orbitrap




AADDTWEPFASGK




AADDTWEPFASGK



YVLTQPPSVSVAPGQTAR



YVLTQPPSVSVAPGQTAR




Interference

Analysis of low abundance proteins is **sensitive to interference** from other components of the sample.



MS1 interference: other components of the sample that overlap with the isotope distribution.

MS/MS interference: other components of the sample with same precursor and fragment masses as the transitions that are monitored.

MS1 interference



Quantitation using MRM







Data taken from CPTAC
Verification Work Group
Study 7.





10 peptides
3 transitions per peptide
Concentrations 1-500 fmol/ μ l
Human plasma background
8 laboratories
4 repeat analysis per lab

Addona et al., Nature
Biotechnol. 27 (2009) 633-641

Quantitation using MRM




Ratios of intensities of transitions



Model: Noise and Interference

Can the knowledge of the relative intensity of the transitions be used to correct for interference?



- Noise is a normally distributed increase or decrease in the intensity.



- Interference is an increase in the intensity of one or more transitions.

Detection of interference

Interference is detected by comparing the ratio of the intensity of pairs of transitions with the expected ratio and finding outliers.

Transition i has interference if $Z_{threshold} < Z_i$

where $Z_{threshold}$ is the interference detection threshold;

$$Z_i = \max_{j \neq i} Z_{ji} = \max_{j \neq i} \frac{r_{ji} - I_j / I_i}{\sigma_{ji}};$$


z_{ji} is the number of standard deviations that the ratio between the intensities of transitions j and i deviate from the noise;

I_i and I_j are the log intensities of transitions i and j ;


r_{ji} is the median of the log intensity of transitions j and i ;


σ_{ji} is the noise in the ratio.




Error in quantitation after correction in presence of noise but no interference



Corrections for interference













Relative noise = 0.2
 Interference in 1 out of 3 transitions
 Relative intensity of transitions: 1:1:1

Error in quantitation after correction in presence of interference and noise



Correction for MS2 interference



Workflow for quantitation with LC-MS

Standardization


Retention time alignment
Mass calibration
Intensity normalization

Quality Control

Detection of problems with
samples and analysis

Quantitation

Peak detection
Background subtraction
Limits for integration in time and mass
Exclusion of interfering peaks



Takeaway Message

- There are many different ways to quantitate proteins - choose the one that is appropriate for your application.
- In general the earlier you can introduce isotopic labels the better the accuracy.
- Always monitor for interference.

Protein quantitation I: Overview (Week 6)
